

Removable Appliance

Instruction Guide



WEARING YOUR REMOVABLE APPLIANCE

Your Removable Appliance (RA) is to be worn **full time**, only being removed to eat, brush your teeth and when playing contact sport or swimming. If it is not worn as directed, the desired result will not be achieved.

INSERTING AND REMOVING YOUR REMOVABLE APPLIANCE

The first few times you fit your RA on your own, it helps to use a mirror to make sure you have it in the correct position.

Once it is in the correct position, push firmly against the plastic to click the appliance into place. When removing your RA, always use one finger on each side to pull down on the metal wires located at the back of the appliance.

SPEECH

When the Removable Appliance is first fitted, it can be challenging to talk and swallow. This will improve quickly as you adapt to the appliance. To speed up the adaptation process you can:

- Take small sips of water from a drink bottle to help with tongue placement during swallowing.
- Read out loud (or simply talk a lot) to help with pronunciation.

WHAT IF I BREAK OR LOSE MY REMOVABLE APPLIANCE?

Please contact us **ASAP** if you misplace or break your RA. Any interruption in your treatment will prolong the time needed to achieve the best results. Please note that replacement of a lost or broken RA will incur an additional fee.




WILL MY TEETH BE SORE?

There is usually some discomfort for the first few days after the RA has been fitted. Mild pain relieving medication can be taken, if required.

CARING FOR YOUR REMOVABLE APPLIANCE

Your RA requires cleaning twice a day when you brush your teeth. A soft toothbrush and toothpaste can be used to brush your appliance. It is particularly important to brush the surfaces of the appliance that sit against your gums or the roof of your mouth. If you feel your RA is in need of a more thorough clean, you may give it an occasional soak in a solution of equal parts white vinegar and cold water.



We will supply you with a case for your RA. On the rare occasion that it is out of your mouth, put it in the case to avoid losing or breaking it. It is a great idea to label your case with your name and phone number.

Always keep your RA away from pets, especially dogs, as they like its taste and smell. They will destroy it!

WHAT CAN I EAT?

While it is generally recommended that your appliance is removed for eating, it is important that a low sugar, well-balanced diet is adhered to.

After eating, and before putting your RA back into your mouth, rinse your mouth with water. This will dislodge any food trapped around your teeth, which will help reduce the risk of developing tooth decay.

IF YOU ARE EXPERIENCING ANY DIFFICULTIES WITH YOUR REMOVABLE APPLIANCE PLEASE CALL OUR CLINIC

Pure Orthodontics

phone: 9370 3155 | www.pureorthodontics.com.au | email info@pureorthodontics.com.au